



Crucibles

STAR® System

Innovative Technology

Versatile Modularity

Swiss Quality

High Quality Crucibles
For Best Thermal Analysis Results

METTLER TOLEDO

Crucibles for Thermal Analysis


The Optimal Type for Every Application

Crucibles serve as containers for samples during DSC and TGA measurements. They guarantee that the sensor is not contaminated during the measurement. The type of crucible used for a measurement can have a large effect on the quality of the results obtained, and also influences important characteristics of the measuring cell. Considering the relevant factors before the measurement can often help save time later on when interpreting the curve.

Features and benefits of the METTLER TOLEDO thermal analysis crucibles:

- **Large volume crucibles** – increase the sensitivity
- **Low mass crucibles** – push the resolution to the limit
- **Highest purity** – no reactions between sample and crucible caused by impure materials
- **Material** – guarantees high thermal conductivity and dimensional stability
- **Flat crucible base** – ensures perfect thermal contact and minimizes artifacts

For crucible selection, the following factors have to be considered:





Choosing the right crucible guarantees the best results for DSC, TGA and TGA/DSC measurements.


► www.mt.com/ta-crucibles

Reliable Automation

24 Hours a Day, Just Like a Swiss Watch

The sample robot is extremely robust and operates reliably 24 hours a day, 7 days a week. All DSC and TGA models can be automated. The sample robot can process up to 34 samples even if every sample requires a different method and a different crucible. The universal gripper can handle all types of METTLER TOLEDO crucibles or even new crucible types that have been defined by the user.

No sample reaction before measurement



The sample robot can remove the protective aluminum lid from the alumina crucibles. It prevents loss or uptake of water of samples in alumina crucibles. It also protects oxygen-sensitive samples from oxidation.

Unique lid-piercing device



The lid piercing device (wasp) is a unique feature. Samples releasing volatiles or being hygroscopic can be enclosed in hermetically sealed crucibles to prevent mass changes prior to analysis. With the lid piercing kit hermetically sealed aluminum crucibles can be opened just before analysis in the DSC or TGA instrument.

Robot also works with 900 μ L crucibles



Even 900 μ L crucibles can be handled by the robot (special grippers with a longer swiveling arm – ArtNo: 30293250 – are necessary). These special grippers are limited to TGA/SDTA (LF/HT) sensors.

For Pt crucible handling, special grippers for Pt (ArtNo: 30479641) are available.



All our DSC and TGA models can be automated. The sample robot can handle any crucible in any sequence.


► www.mt.com/ta-automation

Wide Range of Crucibles

Obtain the Optimal Measurement Results

Crucibles serve as containers for samples during DSC, TGA/DSC and TGA measurements. They prevent the direct contact of the test specimen with the sensor and thus prevent sensor contamination. Furthermore, intentional choice of the crucibles can contribute to optimum measurement results.


Increase crucible volume/sample mass for higher sensitivity



In the above example, the “weak peak” (black curve: about 10 times the noise level of the system) could only be detected because a 100 μ L crucible was used instead of a standard 40 μ L crucible. If larger crucibles are used it is recommended to limit the heating rate to a maximum of 10 K/min as temperature gradients within the sample are expected.


To improve the curve quality, the same experiment was repeated; this time with a reference 100 μ L crucible that was filled with 105 mg buffer solution (red curve).

Improve resolution with low-mass crucibles



Especially in DSC, the crucible impacts the DSC signal response time and thus the temperature or time resolution. The lower the crucible mass, the better the resolution. The low-mass, 20 μ L light aluminum crucible has the best signal time constant, i.e. peaks become sharper and close-lying peaks can be better separated.

Prevent undesired interaction between sample and crucible



Samples can react with the crucible. In most cases this is an unwanted effect. It can be prevented by choosing a suitable inert crucible. For typical reactive samples such as acid or alkaline materials, inert crucibles made of gold, platinum or glass are recommended. Also gold-plated crucibles from aluminum or steel are considered to be inert.

In some cases the interaction between the test specimen and the crucible is of interest (e.g. insulating material/copper).



METTLER TOLEDO offers a wide range of precision designed and manufactured crucibles with different volumes and materials. You can choose which is best suited for your application. Don't forget how important the sample preparation is.

► www.mt.com/ta-sampleprep

Typical Crucibles for DSC

Appropriate for Your Application

Aluminum is the most commonly used crucible material for DSC measurements. The most frequently used crucible type for DSC (and TGA/DSC) measurements is the 40 μ L aluminum (standard) crucible. It is very shallow (low in height) and has a strong flat base, which ensures that temperature gradients are as low as possible by having a well defined thermal contact with the sensor. A variant with pin is also available for the 40 μ L volume, resulting in highest reproducibility (only for manual insertion). All aluminum crucibles can be combined with various lid types, depending on the application.

Aluminum crucibles with various volumes




Aluminum is the most frequently used crucible material for DSC measurements. METTLER TOLEDO offers aluminum crucibles of various volumes from 20 to 160 μ L. The 20 μ L light Al crucible has the shortest signal time constant. However, it cannot be sealed hermetically. The 25 μ L Al crucible is optimized for enthalpy measurements with high reproducibility requirements. It can be sealed hermetically to suppress endothermic evaporation. Larger volume Al crucibles (100 and 160 μ L) are used when large sample masses are needed.

Choice of 3 lid types



The standard lid is used either to hermetically seal the crucible (to suppress evaporation of solvents) or by piercing a hole in the lid to prevent test specimen from spilling over the crucible. In this case, the atmosphere is practically identical to the atmosphere in the furnace. Lids for robot piercing by the piercing device: The lid is pierced before the transfer to the furnace; 3 needle types – 0.1; 0.7; 1 mm diameter – are available. Lids with preperforated 50 μ m-hole: For measurements in a self-generated atmosphere; evaporation of e.g. water is shifted to its boiling point.

Inert crucibles – gold-plated, pure gold or glass



Gold or gold-plated crucibles (40 μ L) are chemically resistant and the optimal choice if strong acids or bases have to be analyzed. These pans are the only non-aluminum types that can be hermetically sealed by cold-welding; maximum pressure: 0.25 MPa (gold) and 0.2 MPa (gold-plated Al).

One type of glass crucible is offered (volume 15 μ L). This crucible cannot be sealed, however, a glass cover disk, 7 mm in diameter, can be applied that closes the pan.

Alternatives to standard 40 µL Al crucibles

Alternative crucible types are described that improve sensitivity, resolution or inertness compared to standard 40 µL Al crucibles.

| Technique | | Crucible selection | | | |
|---------------------------|--|---|--|---|--|
| DSC, HP DSC, TGA/DSC, TGA | | Standard crucible 40 µL aluminum → select the crucible that fits best: | | | |
| Reason Operation mode | | Improvement of sensitivity | Improvement of resolution and reproducibility | Inertness critical | Catalytic effect |
| | | <ul style="list-style-type: none"> 100 µL Al crucible For detection of weak signals. | <ul style="list-style-type: none"> 20 µL Al crucible Maximal separation of close-lying effects (peak separation). Best for flat materials and powders, less suited for liquids. | <ul style="list-style-type: none"> 40 µL Au-plated Al crucible Lid can be hermetically sealed, for reactive samples (e.g. strong acids and strong bases) at temperatures below 350 °C. | <ul style="list-style-type: none"> 40 µL Cu pan Strong catalytic effects, often used for OIT determinations, catalytic effects are desired. |
| Robotic use | | | <ul style="list-style-type: none"> 25 µL Al crucible Lid can be hermetically sealed, maximal separation of close-lying effects (peak separation), for liquids. | <ul style="list-style-type: none"> 40 µL Au crucible Lid can be hermetically sealed, for reactive samples (e.g. strong acids and strong bases) and high temperatures up to 750 °C. | |
| | | | | <ul style="list-style-type: none"> 30 µL or 70 µL Al₂O₃ For TGA/DSC measurements in the high temperature range up to 1600 °C. | |
| | | | | <ul style="list-style-type: none"> 70 µL sapphire crucibles For C_p measurements, for ash testing and wet chemistry applications. | |
| Manual use (DSC only) | | <ul style="list-style-type: none"> 160 µL Al crucible with pin Manual operation, crucible of choice if large quantities must be applied. | <ul style="list-style-type: none"> 40 µL Al crucible with pin Maximal reproducibility. | <ul style="list-style-type: none"> 15 µL glass pan For reactive samples (e.g. strong acids and strong bases). | |

Selection of the aluminum lid types

Alternatives to standard lids are described.

| Technique | | Lid selection | | | |
|---------------------------------------|--|--|---|--|--|
| DSC, HP DSC, TGA/DSC, TGA | | Select lid type to be combined with Al pans (25 µL, 40 µL, 100 µL, 160 µL) | | | |
| Lid type Material | | Standard lid | Lid with 50 µm hole | Piercing lid (sample protection) | No lid |
| | | <ul style="list-style-type: none"> Lid pierced with needle Protects the sensor from material spilling over; atmosphere outside and inside identical. Lid hermetically sealed | <ul style="list-style-type: none"> Lid with 50 µm hole and in combination with standard 40 µL pans Self-generated atmosphere reduces gas exchange, shifts vaporization to the boiling point. | <ul style="list-style-type: none"> Piercing lids in combination with standard 40 µL pans Makes it possible to limit gas exchange to the start of the experiment; no drying or water uptake of the sample during the waiting period. | <ul style="list-style-type: none"> No lid Optimal gas exchange, typically applied for OIT and OOT measurements (particularly for HP DSC). |
| Aluminum, gold-plated lids, gold lids | | <ul style="list-style-type: none"> Hermetically sealed lid Al, Au, gold-plated; no hole, typically for liquids. | | | |

Typical Crucibles for TGA

Accurate TGA Measurements

The most frequently used crucible type for TGA (and TGA/DSC) measurements is the 70 μL alumina (standard) crucible. All alumina crucibles are provided with aluminum oxide lids. Aluminum oxide crucibles can be reused; appropriate cleaning procedures have to be defined (depending on the degree and type of contamination).

Alumina crucibles – from 30 to 900 μL



The most often used crucibles for TGA measurements are the aluminum oxide (or alumina) crucibles. They allow measurements up to 1600 °C. In addition, aluminum oxide is highly chemically inert – side reactions are very seldom. The volume range of our offering from 30 to 900 μL is enormous.

Aluminum lids for alumina crucibles



All alumina crucibles are provided with alumina lids. These lids have a hole. For applications with samples that contain volatile substances, drying out or moisture uptake is critical. For such cases, aluminum lids are available that can be removed by the sample robot before the crucible is transferred to the furnace.


Sapphire crucibles



Sapphire is a very pure form of monocrystalline aluminum oxide. This is the reason why a sapphire crucible is chemically more resistant than an aluminum oxide crucible. The crucible is notably recommended for melting metals such as Fe and Ni.

Alternatives to standard 70 µL Al₂O₃ (alumina) crucibles

Alternative crucible types that improve sensitivity, resolution or inertness are described compared to standard 70 µL Al₂O₃ crucibles.

| Technique | | Crucible selection | | |
|----------------------------------|--|---|--|--|
| TGA, TGA/DSC, DSC | | Standard crucible 70 µL alumina → select the crucible that fits best: | | |
| | |  | | |
| Reason Furnace type | | Improvement of sensitivity | Improvement of resolution and reproducibility | Inertness-critical |
| Small TGA or TGA/DSC furnace | | | • 40 µL (25 µL) Al crucibles Al crucibles for best sensitivity also for TGA up to 640 °C. | • 40 µL Au crucibles Up to 750 °C, Au crucibles if resolution, sensitivity and inertness are critical. |
| | | | • 30 µL Al ₂ O ₃ crucibles Close-lying effects, peak separation. • 30 µL Pt crucibles High temperature range up to 1600 °C. | • 40 µL Au-plated Al crucibles Up to 350 °C. • 70 µL sapphire crucible C _p measurements; measurements of alloys, temperatures up to 1100 °C. |
| TGA: Large or HT furnace | | • 100 µL Al crucibles Al Crucibles for best sensitivity also for TGA up to 640 °C. • 150 µL, 300 µL, 600 µL Al ₂ O ₃ crucibles Weak signal, increase of sample mass required. • 150 µL Pt crucibles Lid not sealed, temperature range > 600 °C, weak signal, reproducibility critical. | • 30 µL Al ₂ O ₃ crucibles Close-lying effects, peak separation. • 30 µL Pt crucibles High temperature range up to 1600 °C. | • 70 µL sapphire crucibles C _p measurements; measurements of alloys, temperatures up to 1600 °C. |
| Large furnace / special grippers | | • 900 µL Al ₂ O ₃ crucibles Special grippers for robot, weak signal, large mass/volume needed (EGA measurements). | | |
| Sapphire disks (TGA) | | | Sapphire disks are always necessary if Pt crucibles are in contact with Pt sensors (DTA). | |

Typical applications intended for Pt crucibles:


- TGA/DSC experiments with a focus on DSC signals at higher temperatures.
- Heat capacity measurements at temperatures > 600 °C.
- Metal characterization measurements (mainly below the melting temperature of the metal sample).
- Humidity measurements with salts.

High-Pressure Crucibles

For the Investigation of Chemical Reactions

METTLER TOLEDO offers a variety of high-pressure crucibles that serve as sample containers for thermoanalytical measurements. These hermetically sealed HP crucibles withstand pressures up to 15 MPa.

High-pressure (HP) crucible – for reuse



The relatively light and flat construction of the HP steel crucible results in low temperature gradients. Compared to larger HP crucibles, the obtained DSC results are better. The thread and the sealing tool with defined torque enable the crucible to be easily and securely sealed. After the measurement, the crucible can be cleaned for reuse (about 20 times, using a new gold-plated copper seal). NiCr (NiCr20TiAl) HP Crucibles: For manual usage and large volumes, two types of NiCr crucibles (270 und 500 μL) are available.

High-pressure crucibles, stainless steel, single-use



These HP gold-plated steel crucibles (25 μL and 40 μL volume) have proven highly effective in the field of safety investigations. These are the optimal single-use crucibles for a pressure up to 15 MPa. The crucible lid is pressed (pressure of about a ton) into the crucible with the aid of a special toggle press so that the rupture disk (gold) hermetically seals the crucible.

Medium-pressure crucibles, single-use



The medium pressure crucible allows a maximum pressure to be generated of about 2 MPa. To use the standard sealing press, the die set (type A3) must be used.

The medium pressure crucibles can be sealed with either a FPM O-ring (that allows water-vapor to permeate) or with a PCTFE O-ring. However, this polymer shows a DSC melting peak at 220 $^{\circ}\text{C}$.

Selection of high-pressure crucibles

The high-pressure portfolio is shown.

| Technique | | Crucible selection | | |
|--|----------|---|---|---|
| DSC, (HP DSC) | | High pressure crucibles | | |
| Material | Lid type | Reusable HP crucibles for moderate reactive samples | Single-use HP crucibles | Reusable, inert HP crucibles |
| | | | | |
| Steel, max. pressure 15 MPa | | <ul style="list-style-type: none"> • 30 µL crucible with burst disc - Suppression of boiling, for highly reactive samples. - Reusable, up to 20 times. | | |
| Gold-plated steel, max. pressure 15 MPa | | | <ul style="list-style-type: none"> • 25 µL, 40 µL crucibles with lid and burst disc - Safety measurements. - Not reusable. | <ul style="list-style-type: none"> • 30 µL crucible with lid ¹⁾ - Suppression of boiling, inert material necessary, highly reactive material. - Reusable, up to 20 times. |
| NiCr ²⁾ alloy, max. pressure 10 MPa | | | | <ul style="list-style-type: none"> • 270 µL, 500 µL crucible with lid ³⁾ and seal - Reusable, no limitation. |

¹⁾ Seals necessary, not included.

²⁾ NiCr: alloy of Ni, Cr, Ti, Al.

³⁾ Furnace expander required; after each measurement, a new seal disk is used.

Crucible material alternatives to aluminum (ambient pressure)

Alternative crucible materials and their advantages compared to the predominantly used aluminum are shown.

| Technique | | Crucible selection | | | |
|----------------------------------|--------|--|--|--|---|
| DSC, HP DSC, TGA/DSC, TGA | | Select crucible material other than aluminum | | | |
| Material | Reason | Platinum (Pt80Rh20) | Al ₂ O ₃ (alumina) | Sapphire | Gold |
| | | | | | |
| Reactivity / chemical resistance | | <ul style="list-style-type: none"> • 30 µL, 70 µL, 150 µL Catalytic effects possible (often undesired). Normally used at temperatures up to 1600 °C. | <ul style="list-style-type: none"> • 30 µL, 70 µL, 150 µL, 300 µL, 600 µL, 900 µL Inert material, easy to handle, Alumina is critical in combination with some metals such as iron. Not optimal for highest resolution and reproducibility. | <ul style="list-style-type: none"> • 70 µL Optimal for C_p measurements and measurements at very high temperatures, sapphire is chemically extremely resistant. Sapphire is the optimal choice for many metals at higher temperatures. | <ul style="list-style-type: none"> • 40 µL Inert material, often used in combination with acid or alkaline samples. |
| Degree of reactivity | | Low and undesired. Pt is the optimal choice for substances that react with Al ₂ O ₃ or sapphire. | Low reactivity (critical with certain metals such as iron). | None (highest degree of inertness). | Low reactivity |

Platinum crucibles: Platinum crucibles are mainly used for TGA/DSC or TGA measurements with end temperatures above 640 °C. SDTA and DSC curves measured with platinum crucibles are usually better than those obtained using crucibles made of alumina, which has a poorer thermal conductivity. In addition, Pt crucibles are often used for accurate C_p determinations (DSC and TGA instruments).

Innovative Accessories for Crucibles

Sealing Press Types and Accessory Box

The standard sealing press allows the pan to be sealed very easily. Under the pressure of the plunger, the pan is cold-welded hermetically with the lid. There is a variety of plungers and dies for various crucible types.

Sealing of many crucible types



A1: Die and plunger set for the aluminum light 20 μL crucibles. Hermetical sealing is not possible with this crucible type.

A2: Die and plunger for all crucibles that can be sealed by cold-welding: 40, 100 and 160 μL aluminum crucibles; 40 μL gold crucibles. All these crucibles can be hermetically sealed.

A1a: Die and plunger for the 40 μL gold-plated aluminum crucibles; they can also be hermetically sealed.

Crucible handling set



The crucible handling set provides a range of tools that are helpful for sample preparation:

- A funnel for filling the pan with sample
- Tweezers for handling of sample, pans and lids
- Standard needle (1 mm) and a thin needle (0.35 mm thick)
- Crucible holder for crucible handling and safe transfer to the instrument.

Sealing of high-pressure crucibles



High-pressure crucibles (25 μL and 40 μL) can be sealed by using a special sealing press (Maeder press). A particular plunger/die set is included for these two types of crucibles.

Die and plunger




A1: Die and plunger for Al crucible light.

ArtNo: 51140547

A1a: Die and plunger for gold-plated aluminum crucibles.


ArtNo: 51142272



20 μL


40 μL

High pressure crucible sealing tool for reusable crucibles



B: This assembly consists of a lower part to keep the hexagonal crucible and an upper part (turning head) that fits on the crucible cover. The crucible is closed with the seal disk by turning until a sliding clutch is activated. After the measurement, the crucible is opened with the same tool.

ArtNo: 51119915




30 μL

Die and plunger



A2: Die and plunger; sealable by cold welding (for hermetical sealing; included in delivery).

ArtNo: 27809 / 27386



25 μL

40 μL

100 μL

160 μL

High pressure crucible sealing tool for reusable NiCr crucibles



C: This is the tool for the NiCr crucibles. They are fixed in one of the tolls and the cover in the other to seal them with the required torque.

ArtNo: 650067



500 μL

270 μL

Die and plunger



A3: Die and plunger for medium pressure crucible.

ArtNo: 119428



120 μL

Tool kit for toggle press



D: This assembly fits in the Maeder press used to seal the 25 and 40 μL high pressure crucibles.

ArtNo: 26733



25 μL

40 μL

Technical Data

| Article no. | Designation | Quantity | Volume [μL] | With pin | Without pin | With lid | Without lid | Pan weight [mg] | Lid weight [mg] | Max. pressure [MPa] | Max. temperature [°C] | Ø Base [mm] | Height [mm] without lid |
|--------------|---|------------|-------------|----------|-------------|----------|-------------|-----------------|-----------------|---------------------|-----------------------|-------------|-------------------------|
| 26763 | Aluminum crucible (standard) | 100 | 40 | | • | • | | 33 | 17 | 0.2 | 640 | 6 | 1.6 |
| 51119810 | Aluminum crucible light | 100 | 20 | | • | • | | 14 | 9 | 0.2 | 640 | 6 | 1.6 |
| 30085850 | Aluminum crucible | 100 | 25 | | • | • | | 28 | 17 | 0.2 | 640 | 4 | 1.6 |
| 27331 | Aluminum crucible standard | 100 | 40 | • | | • | | 33 | 17 | 0.2 | 640 | 6 | 1.6 |
| 51119870 | Aluminum crucible standard | 400 | 40 | | • | | • | 33 | – | 0.2 | 640 | 6 | 1.6 |
| 51141448 | Aluminum crucible / light lid | 100 | 40 | • | | • | | 33 | 9 | 0.2 | 640 | 6 | 1.6 |
| 51143092 | Aluminum crucible large | 100 | 100 | | • | | • | 63 | – | 0.2 | 640 | 6 | 4.2 |
| 51119872 | Aluminum crucible large | 400 | 100 | | • | | • | 63 | – | 0.2 | 640 | 6 | 4.2 |
| 27811 | Aluminum crucible large | 40 | 160 | • | | • | | 80 | 17 | 0.2 | 640 | 6 | 6.4 |
| 51119871 | Aluminum lid standard | 400 | lid | | | | | – | 17 | | 640 | | |
| 51140832 | Aluminum lid pierced 50 μm | 400 | lid | | | | | – | 17 | | 640 | | |
| 51119873 | Aluminum piercing lid | 400 | lid | | | | | – | 16 | | 640 | | |
| 51142973 | Aluminum crucible gold plated | 10 | 40 | | • | • | | 70 | 40 | 0.2 | 350 | 6 | 1.6 |
| 27220 | Gold crucible | 6 | 40 | | • | • | | 270 | 130 | 0.25 | 750 | 6 | 1.65 |
| 51140407 | Copper crucible | 100 | 40 | | • | | • | 75 | – | | 750 | 6 | 1.65 |
| 17780 | Glass crucible | 5 | 15 | | | | | 55 | | | 500 | 7 | 1 |
| 17782 | Glass cover disk (7 mm) | 5 | disk | | | | | | 15 | | | 7 | 0.2 |
| 24123 | Alumina crucible (standard) | 20 | 70 | | • | • | | 185 | 60 | | 2000 | 6 | 4.5 |
| 51140843 | Alumina crucible small | 20 | 30 | | • | • | | 95 | 60 | | 2000 | 6 | 2.6 |
| 24124 | Alumina crucible large | 20 | 150 | | • | • | | 280 | 110 | | 2000 | 8 | 4.5 |
| 30267108 | Alumina crucible large | 10 | 300 | | • | • | | 520 | 100 | | 2000 | 8 | 9 |
| 30077260 | Alumina crucible large | 4 | 600 | | • | • | | 880 | 230 | | 2000 | 12 | 9.2 |
| 51119960 | Alumina crucible large | 4 | 900 | | • | • | | 770 | 270 | | 2000 | 12 | 10 |
| 51119649 | Aluminum lid for Al ₂ O ₃ 30/70 | 40 | lid | | | • | | | 50 | | | | |
| 51140477 | Aluminum lid for Al ₂ O ₃ 150/300 | 40 | lid | | | • | | | 75 | | | | |
| 30077266 | Aluminum lid for Al ₂ O ₃ 600 | 40 | lid | | | • | | | 205 | | | | |
| 51140469 | Aluminum lid for Al ₂ O ₃ 900 | 40 | lid | | | • | | | 185 | | | | |
| 51140842 | Platinum crucible small | 4 | 30 | | • | • | | 180 | 85 | | 1600 | 6 | 2.3 |
| 51119654 | Platinum crucible medium | 4 | 70 | | • | • | | 275 | 85 | | 1600 | 6 | 4.2 |
| 24126 | Platinum crucible large | 4 | 150 | | • | • | | 420 | 140 | | 1600 | 7.3 | 4.2 |
| 51140845 | Sapphire crucible | 4 | 70 | | • | • | | 200 | 60 | | 2000 | 6 | 4.5 |
| 51140404 | High pressure crucible | 3 | 30 | | • | • | | 300 | 270 | 15 | 750 | 7 | 2.5 |
| 51140405 | High pressure crucible gold plated | 3 | 30 | | • | • | | 325 | 270 | 15 | 350 | 7 | 2.5 |
| 51140403 | Seal for high pressure crucible | 60 | seal | | | | | – | – | | | | |
| 30077139 | High pressure crucible small | 25 | 25 | | • | • | | 620 | 280 | 15 | 400 | 6.4 | 4.7 |
| 26732 | High pressure crucible / pin | 25 | 40 | • | | • | | 1150 | 280 | 15 | 400 | 7 | 5.9 |
| 26731 | High pressure crucible | 25 | 40 | | • | • | | 1150 | 280 | 15 | 400 | 7 | 5.9 |
| 650072 | High pressure crucible | 1 | 270 | • | | • | | 1060 | 1070 | 10 | 750* | 7.6 | 8.4 |
| 650066 | High pressure crucible | 1 | 500 | • | | • | | – | – | 10 | 750* | 7.6 | 14.4 |
| 27216 | Seal for high pressure crucible | 1 | seal | | | | | – | – | | 750* | | |
| 26929 | Medium pressure crucible | 25 | 120 | • | | • | | 170 | 140 | 2 | 250 | 6 | 5.5 |
| 29990 | Medium pressure crucible | 25 | 120 | | • | • | | 170 | 140 | 2 | 250 | 6 | 5.5 |
| 26933 | Seal for medium pressure crucible | 30 | seal | | | | | – | – | | 230 | | |

* Above 500 °C and in the presence of oxygen, scaling of the material may occur.

| Article no. | Designation | Material | Sealing tool | For DSC | For TGA | Suitable for sample robot | Turn table U1 | Turn table U2 | Turn table U3 |
|--------------|---|--|--------------|---------|---------|---------------------------|---------------|---------------|---------------|
| 26763 | Aluminum crucible (standard) | Al 99.99% | A2 | • | • | • | | • | • |
| 51119810 | Aluminum crucible light | Al 99.99% | A1 | • | • | • | • | • | • |
| 30085850 | Aluminum crucible | Al 99.99% | A2 | • | • | • | • | • | • |
| 27331 | Aluminum crucible standard | Al 99.99% | A2 | • | | | | | |
| 51119870 | Aluminum crucible standard | Al 99.99% | A2 | • | • | • | • | • | • |
| 51141448 | Aluminum crucible / light lid | Al 99.99% | A2/A1 | • | | | | | |
| 51143092 | Aluminum crucible large | Al 99.99% | A2 | • | • | • | • | • | • |
| 51119872 | Aluminum crucible large | Al 99.99% | A2 | • | • | • | • | • | • |
| 27811 | Aluminum crucible large | Al 99.99% | A2 | H | | | | | |
| 51119871 | Aluminum lid standard | Al 99.99% | | • | • | • | • | | • |
| 51140832 | Aluminum lid pierced 50 µm | Al 99.99% | | • | • | • | • | | • |
| 51119873 | Aluminum piercing lid | Al 99.99% | | • | • | • | • | | • |
| 51142973 | Aluminum crucible gold plated | Al 99.5%, 5 µm gold plated | A1a | • | • | • | • | • | • |
| 27220 | Gold crucible | Au 99.99% | A2 | • | • | • | • | • | • |
| 51140407 | Copper crucible | E-Cu 99.90% | | • | • | • | • | • | • |
| 17780 | Glass crucible | Borosilicate glass | | • | • | | | | |
| 17782 | Glass cover disk (7 mm) | Borosilicate glass | | • | • | | | | |
| 24123 | Alumina crucible (standard) | Al₂O₃ 99.7% | | • | • | • | • | • | • |
| 51140843 | Alumina crucible small | Al ₂ O ₃ 99.7% | | • | • | • | • | • | • |
| 24124 | Alumina crucible large | Al ₂ O ₃ 99.7% | | | E | • | • | | • |
| 30267108 | Alumina crucible large | Al ₂ O ₃ 99.7% | | | E | • | • | • | • |
| 30077260 | Alumina crucible large | Al ₂ O ₃ 99.5% | | | E | • | • | • | • |
| 51119960 | Alumina crucible large | Al ₂ O ₃ 99.8% | | | E | G | | | |
| 51119649 | Aluminum lid for Al ₂ O ₃ 30/70 | Al 99.5% | | • | • | | | | |
| 51140477 | Aluminum lid for Al ₂ O ₃ 150/300 | Al 99.5% | | | E | | | | |
| 30077266 | Aluminum lid for Al ₂ O ₃ 600 | Al 99.5% | | | E | | | | |
| 51140469 | Aluminum lid for Al ₂ O ₃ 900 | Al 99.5% | | | E | | | | |
| 51140842 | Platinum crucible small | Pt Rh 20% | | • | • | • | • | • | • |
| 51119654 | Platinum crucible medium | Pt Rh 20% | | • | • | • | • | • | • |
| 24126 | Platinum crucible large | Pt Rh 20% | | • | E | • | • | | • |
| 51140845 | Sapphire crucible | α-Al ₂ O ₃ 99.9% | | • | • | • | • | • | • |
| 51140404 | High pressure crucible | X2 CrNiMo18143 | B | • | | • | | | • |
| 51140405 | High pressure crucible gold plated | X2 CrNiMo18143, 5 µm gold plated | B | • | | • | | | • |
| 51140403 | Seal for high pressure crucible | Cu, 2 µm gold plated | | | | | | | |
| 30077139 | High pressure crucible small | X2 CrNiMo18143, 5 µm gold plated | D | • | | • | | | |
| 26732 | High pressure crucible/pin | X2 CrNiMo18143, 5 µm gold plated | D | • | | | | | • |
| 26731 | High pressure crucible | X2 CrNiMo18143, 5 µm gold plated | D | • | | • | | | • |
| 650072 | High pressure crucible | NiCr 20 TiAl | C | F | | | | | |
| 650066 | High pressure crucible | NiCr 20 TiAl | C | F | | | | | |
| 27216 | Seal for high pressure crucible | Au 700/531 | | | | | | | |
| 26929 | Medium pressure crucible | X5 CrNi 18 9 | A3 | • | | | | | |
| 29990 | Medium pressure crucible | X5 CrNi 18 9 | A3 | • | | • | | • | • |
| 26933 | Seal for medium pressure crucible | PCTFE | | | | | | | |

E: Only with large furnace (900/600 µL alumina crucible: only without the lid piercing kit)

F: Only with DSC20, 25, 27HP, 30; STAR® DSC82x and STAR® System DSC 1/2/3 with furnace expander (without sample changer)


G: Only with the special gripper (30293250) for sample robot

H: The crucible can be used without furnace expander if the lid is reversed

Extremely Wide Application Range Crucibles for Every Sample Type

Good thermal analysis measurements also depend also on the correct choice of crucibles and lids. Lids seal the crucible hermetically to avoid evaporation of the sample and to avoid interferences with the surrounding atmosphere. In contrast, open crucibles without a lid (or with a lid with a big hole) allow the ambient atmosphere to come into contact with the sample.

Sealing possibilities






Crucibles are fundamental for the studies of materials such as thermoplastics, thermosets, elastomers, composites, adhesives, foodstuffs, pharmaceuticals and chemicals.


► www.mt.com/ta-applications

Influence of crucible sealing





This example describes the influence of crucible sealing on the resulting measurement curves obtained from a DSC experiment. Two polyamide samples were measured using standard 40 μ L aluminum crucibles; one was left open (see the green curve), and the other was hermetically sealed (see the blue curve). At room temperature, the moisture in the sample begins to evaporate in the open crucible; this results in a broad endothermic peak on the green DSC curve and the glass transition at about 30 $^{\circ}$ C is hardly detectable. In contrast, evaporation is suppressed in the hermetically sealed crucible and the glass transition can be easily identified.

Cholesteryl myristate measured with various crucibles




This example shows the DSC curves of cholesterol myristate – a liquid crystal that undergoes two liquid-liquid transitions after melting. The only difference between the measurements was the type of crucible used. The blue curve measured using a 30 μ L gold plated high-pressure crucible shows a broad melting peak at about 74 $^{\circ}$ C, followed by two weak effects corresponding to a liquid-liquid transition at about 81 and 86 $^{\circ}$ C, respectively. These latter two effects exhibit sharper peaks when the sample is measured in a 20 or 40 μ L aluminum crucible.

Sapphire crucibles used to study iron




In this example, the phase transitions of iron are illustrated. The measurements were performed in argon using a $70 \mu\text{L}$ sapphire crucible. These crucibles are of Al_2O_3 material that has been crystallized at very high temperature and is resistant towards molten iron (in contrast to conventional alumina crucibles). The examples shows 3 modifications of iron (α -modification prevalent below 920°C ; γ -modification formed at about 920°C and a δ -modification formed at about 1397°C). The iron melts at about 1540°C .

Inert crucibles suppress undesired interactions




For most applications, a reaction between the crucible material and the sample is undesired. Aluminum crucibles are suitable for most samples but may react with certain acids or alkaline materials. In this example, a wood pulp specimen containing sodium hydroxide and sodium sulfide reacted with the aluminum crucible to produce a large exothermic peak. This is clearly illustrated by the green curve. In contrast, when the same sample is measured in an inert crucible such as alumina or gold, the desired evaporation peak of the water is clearly observed as an endothermic peak; see red and blue curves.

Copper crucibles enhance catalytic effect




OIT experiments – short for oxidative induction time – provide information about the oxidative stability of a sample. Here, a high-density PE sample was heated to 200 °C in an inert atmosphere. After 12.5 min at 200 °C, the gas is switched from nitrogen to oxygen. The time from the gas switching to the onset of the oxidation reaction is defined as OIT. The results show a big difference in the stability of the sample with respect to the pan material. The sample with the copper pan has a significantly lower OIT than the sample with the aluminum pan as copper increases the reaction time. Such experiments are important for the cable industry (ASTM D3895).

Platinum crucibles for c_p measurements




Using a TGA/DSC, this example illustrates the determination of the specific heat capacities (C_p) of two metals (molybdenum and nickel; blue line) and the comparison with the literature value (black line with symbol). For C_p measurements, crucibles of platinum are recommended as they have an excellent thermal conductivity. In general, the crucible should be used with lid and its size should be as small as possible. As Pt can stick to Pt sensors at temperatures around 1100 °C and above, we recommend placing sapphire disks between the crucibles and sensor.

Gypsum measured with different pan-lid combinations



The type of sealing also influences a DSC measurement. This can be illustrated using calcium sulfate dihydrate, commonly referred to as gypsum. When gypsum is heated, dehydration occurs in two steps to form hemihydrate followed by anhydrite. These dehydration steps were studied in an open crucible and two sealed crucibles, one pierced with a 1 mm hole and the other with a 50 µm hole. As shown here, different types of sealings have a significant influence on the results. The sample measured in an open crucible (top curve) exhibits one broad endothermic peak due to dehydration. When the same sample is measured in sealed crucibles with pierced lids, the two dehydration steps become visible (see middle and bottom curves). The different phases are best separated when a very small hole is pierced in the lid – in this case 50 µm.

Moisture determination in a polymer membrane



When a low moisture polymer membrane sample is heated in a TGA, a large amount of sample is needed to accurately determine its moisture content. Therefore, about 400 mg of sample was measured in a 900 µL crucible. Between sample preparation and the actual measurement, the sample loses about 0.05% of its initial mass due to evaporation. Upon heating to 130 °C, an additional 0.1% is lost. This corresponds to about 400 µg of sample. In contrast, a 70 µL crucible could only support an initial sample mass of about 30 mg, and the overall mass loss step would correspond to 45 µg. This results in a measurement uncertainty of about 10% in contrast to only 0.3% when using a 900 µL crucible filled with 400 mg of sample.

Simple, Intuitive Operation

Straightforward, Efficient and Secure

STAR[®] is the most comprehensive thermal analysis software on the market and provides unrivalled flexibility and unlimited evaluation possibilities. The flexible STAR[®] software consists of the base software and a large number of application-specific options to satisfy your current and future requirements.

Crucible definition in Install Window


Crucible types are defined in the Install Window and belong to every method. The STAR[®] software allows users to define their own crucible types with new dimensions.

FlexCal™: Automatically get the correct adjustment parameters



- FlexCal database for the highest accuracy: STAR[®] software stores the complete adjustment data in the database for every crucible, gas and module combination.
- FlexCal factors for fast results: For the case that there is no specific adjustment data available, the system takes a base adjustment dataset and calculates a dataset for another crucible, gas and module combination.


TabVal calibration and adjustment



Calibration and adjustment has been revised and optimized so that the possibilities offered by FlexCal are more advantageous.

For each calibration substance, there is now a result value with details of the tolerance for the particular heating rate. The evaluation of the adjustment or calibration measurement is fully automatic.

Complete Thermal Analysis System



A complete thermal analysis system consists of the basic six complementary measuring techniques, each of which brings fast and accurate results. Additional knowledge can be obtained by means of several hyphenated techniques.

► www.mt.com/ta-software

World-Class Service and Support Provide Results You Can Trust

METTLER TOLEDO's portfolio of services is designed to ensure the continuous performance and reliability of your thermal analysis systems. Factory-trained in Switzerland, our worldwide teams bring the professional expertise and know-how needed to provide you with the highest level of after-sales support, as well as the experience necessary to optimize services for your own particular needs.

Large library of overview and how-to videos



Increase your productivity by learning from our video libraries and more than 600 applications.

Videos

► www.mt.com/ta-videos

Handbooks

► www.mt.com/ta-handbooks

Applications

► www.mt.com/ta-applications

Comprehensive training courses



We offer effective classroom training sessions.

► www.mt.com/ta-training

For self-training purposes, you may purchase the tutorial kit consisting of 23 well-chosen application examples and their corresponding test substances.

► www.mt.com/ta-tutorial

Register for an e-training course

► www.mt.com/ta-ettraining

Bi-annual application magazine



Every year, thermal analysis generates a large number of scientific results and discoveries. Interesting examples from different application fields and industries are published in our UserCom magazine.

► www.mt.com/ta-usercoms

Thermal Analysis eNewsletter

Sign up to receive quarterly news on the latest trends in thermal analysis, including applications, webinars, courses and how-to-videos.

► www.mt.com/ta-knowledge



The EduPac consists of detailed checklists that provide the sales persons tools to define the training content necessary for a certain customer. In addition, the EduPac helps to estimate the duration of such a customer training.

► www.mt.com/ta-GTAP

For More Information

Outstanding Services

METTLER TOLEDO offers you valuable support and services to keep you informed about new developments and help you expand your knowledge and expertise.

► www.mt.com/ta-services

News on Thermal Analysis

Informs you about new products, applications and events.

► www.mt.com/ta-news

► www.mt.com/ta-app

► www.mt.com/ta-knowledge

UserCom

Our popular, biannual technical customer magazine, where users and specialists publish applications from different fields.

► www.mt.com/ta-usercoms

Applications

If you have a specific application question, you may find the answer in the application database.

► www.mt.com/ta-applications

Webinars

We offer web-based seminars (webinars) on different topics.

► www.mt.com/ta-webinars (Live Webinars)

► www.mt.com/ta-ondemand (On Demand Webinars)

Handbooks

Written for thermal analysis users with background information, theory and practice, useful tables of material properties and many interesting applications.

► www.mt.com/ta-handbooks

Videos

Our technical videos explain complex issues concerning thermal analysis instrumentation and the STAR[®] software – whether it's sample preparation, installation, creating experiments or evaluating measurement results.

► www.mt.com/ta-videos

Training

Classroom training is still one of the most effective ways to learn. Our User Training Courses will help you get the most out of your equipment. We offer a variety of one-day theory and hands-on courses aimed at familiarizing you with our thermal analysis systems and their applications.

► www.mt.com/ta-training (Classroom)

► www.mt.com/ta-ettraining (Web-based)

www.mt.com/ta

For more information

METTLER TOLEDO Group

Analytical Instruments

Local contact: www.mt.com/contacts

Subject to technical changes

© 12/2022 METTLER TOLEDO

All rights reserved. 30208680H

Marketing MatChar / MarCom Analytical



Quality certificate. Development, production and testing according to ISO 9001.



Environmental management system according to ISO 14001.



"European conformity". The CE conformity mark provides you with the assurance that our products comply with the EU directives.